



Key Features

- Alternative for LN2 Storage
- Convenience
- Reliable Performance

More Detail:

http://www.haiermedical.com/products/preservation/150_cryo_freezer/dw-150w200.html

Specifications

Technical Data

Type	-150°C Cryo Freezer
Cabinet Type	Chest
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	HFC
Noise Level (dB)	56

Performance

Cooling performance(°C)	-150
Temperature Range(°C)	-126~-150

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	380/50
--------------------	--------

Power(W)	2800
Electrical Current(A)	11
Dimensions	
Capacity(L/Cu.Ft)	200/7.0
Net/Gross Weight(approx)	315/385(kg) 694.4/848.8(lbs)
Interior Dimensions(W*D*H)	667*462*650(mm) 26.3*18.2*25.6(in)
Exterior Dimensions(W*D*H)	1670*850*1060(mm) 65.7*33.5*41.7(in)
Packing Dimensions(W*D*H)	1810*940*1240(mm) 71.3*37.0*48.8(in)
Container load (20'/40'/40'H)	6/13/13
Functions	
High/Low Temperature	Y
Remote Alarm	Y
Hot Condenser	Y
Power Failure	Y
High/Low Voltage	Y
Sensor Error	Y
Low Battery	Y
High Ambient Temp	Y
Door Ajar	N/A
Accessories	
Caster	Y
Foot	Y
Test Hole	Y
Shelves/Inner Doors	-/2
Temperature Recorder	Optional
USB Interface	N/A
Rs232/485 Interface	Optional
CO2 Backup System	Optional
LN2 Backup System	Optional
Racks And Boxes	Optional

CO2 Backup System

LN2 Backup System

Chart Recorder

Others

Certificate
