



INTRODUCING THE NEW MAGNETIC STIRRER PLATES

SAMCO magnetic stirrers mark a new generation in stirrer plate design. Ultra thin, at just 13.2mm and 13.4mm deep, the X1-150 and X4-450 models work with a breakthrough motorless stirrer mechanism, so there are no moving parts to require maintenance.

Product features include:

- Smart, slim-line, compact design
- Digital display, with speed selection and timer functions
- Speed selection from 15 to 1,200 rpm (1,500 rpm for model X1-150)
- Continuous mode & timer settings from 1-99 minutes
- Pulse mode option - reverses the motion direction every 30 minutes for thorough mixing
- Micro controller programming software enables last run memory and allows gradual stirring acceleration
- Stainless steel, corrosion resistant working area
- IP 65 compliant design (water spillage resistant)

Electrical features include:

- 24V, low voltage device for safety
- Universal power supply (95-230V - 50/60Hz), with interchangeable pins supplied

Model X1-150 Stirrer Plate

Single station

Plate/platform size: 102.5mm x 107mm



Model X1-150



SELLING PRICE:
RM 600.00 /unit

Code	Dimensions mm	Stirring capacity ml	Speed setting rpm
M900/01	156 x 110 x 13.2	800	15 to 1,500

