



### Features

- ▶ Heating element is coated by glass fiber, flexible and easy to coil
- ▶ As flexible heating tapes and cords, easy to use for glass tube, column, pipelines, etc
- ▶ Heating cords may be most suitable when material is small size and complicated.
- ▶ Max. temperature 450°C

### Specifications

Model	WidthXLength (mm)	Watts	Price	Model	WidthXLength (mm)	Watts	Price	Model	WidthXLength (mm)	Watts	Price
HT 2506	25×600	150		HC 1206	12×600	80		HC 0706	∅7×600	80	
HT 2510	25×1000	250		HC 1210	12×1000	120		HC 0710	∅7×1000	120	
HT 2515	25×1500	380		HC 1215	12×1500	180		HC 0715	∅7×1500	180	
HT 2520	25×2000	500		HC 1220	12×2000	240		HC 0720	∅7×2000	240	
HT 2525	25×2500	620		HC 1225	12×2500	300		HC 0725	∅7×2500	300	
HT 2530	25×3000	750		HC 1230	12×3000	360		HC 0730	∅7×3000	360	

Note : Temperature controller is not included  
 (See page 98~101 for temperature controllers)  
 Order-made: Non-standard goods

How to Order			
1	Size	4	Working condition
2	Capacity	5	Electric capacity
3	Working temperature	6	Power supply

\* Please give us above mentioned information

### Features

- ▶ Heating tapes for high temperature : Max. 600°C

### Specifications

Model	WidthXLength (mm)	Watts	Price
ST 2506	25 × 600	250	
ST 2510	25 × 1000	380	
ST 2515	25 × 1500	560	
ST 2520	25 × 2000	750	
ST 2525	25 × 2500	950	
ST 2530	25 × 3000	1200	

Note : Temperature controller is not included  
 (See page 98~101 for temperature controllers)  
 Order-made on the request of non-standard product