


MULTIPLEX SHOULDN'T BE COMPLEX


Azure Imaging Systems

The Next Generation

Designed for easy Total Protein Normalization, Azure Imaging Systems bring even greater performance and streamlined user interface to the next evolution of quantitative Western blot imaging.



Simultaneous detection of total protein content with AzureRed and two-color NIR Western Blot using AzureSpectra secondary antibodies.




Digital image of 3-color Western Blot using Azure Biosystems 600 imager.




Big performance, small footprint, incomparable value.

Great science starts with high-quality data, and when it comes to imaging gels, blots, and plates, high quality data starts with the Azure Imaging Systems.


Chemi with Azure Red for TPN



4 Channel Fluorescent Western Blot Imaging



User Interface Designed for Quantitative Westerns



Azure Imaging Systems

A complete, flexible imaging system line, designed to provide the flexibility needed for your applications.

Resolution and Applications	200	280	300	400	500	600
Resolution	5.4 MP	6 MP	9 MP	9 MP	9 MP	9 MP
UV Gel Imaging	✓	✓	✓	✓	✓	✓
Trans White Imaging	✓	✓	✓	✓	✓	✓
EPI Visible Color Imaging	✓	✓	✓	✓	✓	✓
Chemiluminescent Westerns		✓	✓	✓	✓	✓
Multiplex Visible Fluorescent Westerns				✓		✓
Multiplex NIR Fluorescent Westerns					✓	✓
Channel for TPN			Optional	✓	Optional	✓



www.azurebiosystems.com • info@azurebiosystems.com

Copyright © 2021 Azure Biosystems. All rights reserved. The Azure Biosystems logo and Azure Biosystems® are trademarks of the Company. All other trademarks, service marks and trade names appearing in this brochure are the property of their respective owners.

MM-0041 R1